




StarMed Ventumask 30

For CPAP therapy with Venturi flow driver, adjustable PEEP valve and pressure manometer



Critical Care ▪ Patient Interfaces

StarMed Ventumask 30 with Venturi flow driver and adjustable PEEP valve and pressure manometer

The Ventumask 30 is a unique design for CPAP therapy that is simply activated by connecting to an appropriate oxygen supply. It can be used in a hospital environment or pre-hospital use, for example in an ambulance.

Advantages and Technical Features

- Ventumask 30 is a complete and easy to use system, where the oxygen supply tubes, a Venturi flow driver and a PEEP valve are integral to a CPAP mask fitted with a pressure gauge
- The built-in adjustable PEEP valve can provide patient pressures in a range of settings up to 20 cmH₂O
- The integrated Venturi flow driver can generate patient flows up to 80 L/min while consuming a reduced amount of compressed oxygen
- O₂% can be adjusted from 30% to 100% at each PEEP level
- The user is guided by a performance table integral to the product, for independent PEEP and O₂% settings

The kit is composed of:
 Venturi flow driver integral to the mask for flows up to 80 L/min

CPAP mask with safety valve and headgear

Integrated pressure gauge

Adjustable PEEP valve for settings up to 20 cmH₂O

Integral oxygen tubes

Dual flow meter Code: STS.15/30
 Dual oxygen flow meter 15 + 30 L/min with O₂ coupling

Ventumask 30 for CPAP therapy with Venturi flow driver, adjustable PEEP valve and pressure manometer

Requirement		Setting	Flow delivered
O ₂ % approx	PEEP cmH ₂ O	Flow O ₂ L/min A + B	Flow approx L/min
30%	5	12 + 0	78
	7.5	12 + 0	75
	10	12 + 0	70
	12.5	12 + 0	65
	15	12 + 0	60
	20	14 + 0	60
40%	5	11 + 8	78
	7.5	11 + 8	72
	10	12 + 6	70
	12.5	12 + 4	65
	15	12 + 4	60
	20	14 + 2	60

Requirement		Setting	Flow delivered
O ₂ % approx	PEEP cmH ₂ O	Flow O ₂ L/min A + B	Flow approx L/min
50%	5	10 + 18	80
	7.5	10 + 18	77
	10	10 + 16	72
	12.5	10 + 14	65
	15	10 + 12	60
	20	12 + 10	60
60%	5	9 + 26	78
	7.5	9 + 26	75
	10	10 + 24	70
	12.5	10 + 22	65
	15	10 + 22	60
	20	10 + 18	60
100%	Feed with a minimum O ₂ flow of 40 L/min		

Code	Description	Size	Box Qty.
CV204S	Ventumask 30 CPAP mask with integral Venturi flow driver and adjustable PEEP valve	Small	5
CV204M/L	Ventumask 30 CPAP mask with integral Venturi flow driver and adjustable PEEP valve	Medium / Large	5
CV204XL	Ventumask 30 CPAP mask with integral Venturi flow driver and adjustable PEEP valve	Extra Large	5
STS.15/30	Dual oxygen flow meter 15 + 30 L/min with O ₂ coupling		1



www.intersurgical.com.au/info/starmed

Make an enquiry

IS1.43_AU • Issue 2 06.16



Unit 4, 151 Beauchamp Road, Matraville, NSW, 2036

T: +61 (0) 2 8048 3300 info@intersurgical.com.au www.intersurgical.com.au

UK ▪ Ireland ▪ France ▪ Germany ▪ Spain ▪ Portugal ▪ Italy
 Netherlands ▪ Sweden ▪ Denmark ▪ Lithuania ▪ Russia
 Czech Republic ▪ Turkey ▪ South Africa ▪ China ▪ Japan
 Taiwan ▪ Philippines ▪ USA ▪ Colombia ▪ Australia



Intersurgical Ltd is certified to ISO 14001, ISO 9001 and ISO 13485

Please think before you print
 Save energy and paper.
 If you must print this e-mail
 please print it double sided