



ClearLite™ anaesthetic face masks

A clear vision for the future 



Anaesthesia ▪ Anaesthetic Face Masks



ClearLite™ anaesthetic face masks

A range of anaesthetic face masks which includes seven sizes to meet the clinical needs of all patient types and sizes. The range is available in sizes 0 neonatal to size 6 large adult.

Due to the enhanced anatomical design, the ClearLite mask offers superior seal performance with minimum pressure, reducing the risk of unnecessary discomfort or potential

trauma. The superior seal will also assist the clinician in providing optimum manual ventilation of the patient in both the anaesthetic and emergency situation.

ClearLite is non-PVC and phthalate free which significantly minimises the impact upon the environment.

Lightweight mask

Provides enhanced grip and comfort in use for the clinician

Anatomical shape

Provides a precise fit for different face types and sizes. Minimum pressure required to achieve an effective seal reducing the risk of inadvertent trauma and nerve damage

Crystal clear body

Provides excellent visibility for confirmation that the patient is well perfused

Flexible non-inflatable cuff

Provides a comfortable and effective airtight seal



Safe and secure connection

This is achieved between the mask and the breathing system in use. The ISO 22mm tapered connection complies with ISO 5356

Reduced mask volume

Reduces potential rebreathing and risks associated with elevated levels of end tidal CO₂

Colour-coding

The comfortable, flexible seal is colour coded to allow quick and easy identification of the correct sized mask

Reduce waste

Available without a hook ring

Video available at www.intersurgical.com.au

Low environmental impact

Non-PVC and phthalate free



Code	Code without hook ring	Description	Box Qty.
7296000	7296001	ClearLite, anaesthetic face mask, size 6, extra large adult, red seal, 22F	20
7295000	7295001	ClearLite, anaesthetic face mask, size 5, large adult, orange seal, 22F	35
7294000	7294001	ClearLite, anaesthetic face mask, size 4, adult, green seal, 22F	35
7293000	7293001	ClearLite, anaesthetic face mask, size 3, small adult, yellow seal, 22F	35
7292000	7292001	ClearLite, anaesthetic face mask, size 2, paediatric, white seal, 22F	25
-	7291000	ClearLite, anaesthetic face mask, size 1, infant, grey seal, no hook ring, 15M	30
-	7290000	ClearLite, anaesthetic face mask, size 0, neonate, blue seal, no hook ring, 15M	40

Make an enquiry



Intersurgical Australia Pty Ltd, 4/151 Beauchamp Road, Matraville, NSW, 2036, Australia

T: +61 (0)2 8048 3300 F: +61 (0)2 8079 0769 info@intersurgical.com.au www.intersurgical.com.au



The manufacturer Intersurgical Ltd is certified to ISO 9001:2015, ISO 13485:2016 and ISO 14001:2015

Please think before you print
Save energy and paper.
If you must print this information sheet please print it double sided.

IS5.31 AU • Issue 3 02.20