




# Anaesthetic Modular Breathing Systems

Prioritising patient safety, providing clinical choice



**Anaesthesia** ▪ Anaesthetic Breathing Systems

## A new standard of care

Intersurgical's Anaesthetic Modular Breathing Systems are a new range of standard breathing systems designed to help improve safety and meet the majority of clinical situations whilst improving inventory management and reduces waste.

All systems benefit from an integrated monitoring line and Silver Knight™ anti-microbial additive improving patient safety and reducing the risk of cross contamination in the operating room.

A choice of either Flextube™ or UniFlow™ breathing systems in a range of tube lengths from 1.6m to 5.6m and the option of Compact™ or Flextube™ bagging limbs. The Anaesthetic Modular Breathing Systems provide maximum product flexibility whilst minimising unnecessary waste and providing improved security of supply.

For ease of reference the new Anaesthetic Modular Breathing Systems uses the last 2 digits of each product code to indicate the length of the breathing system or limb. For example, 2919016 is 1.6m in length


### A Select the type and length of breathing system required

#### Uniflow™ Silver Knight™ breathing systems – available in 4 lengths

[▶ Watch the video](#)

##### Integral monitoring line

Minimises risk of disconnection or the possibility of the monitoring line becoming caught/entangled



##### Large 30mm diameter outer lumen

Ensuring that the resistance to flow of the expired gases is kept to a minimum

##### Each UniFlow system comes complete with a leak tester

Easy to test the integrity of the inner lumen as recommended by the AAGBI<sup>1</sup>

##### Push and Twist connectors

Ensures a secure, leak free fit every time



##### Silver Knight

Anti-microbial protection, reducing bacteria by 99.9%





## Flextube™ Silver Knight™ breathing systems – available in 7 lengths


▶ Watch the video

### Reduced connections

Minimises the risk of disconnection or possibility of the monitoring line becoming caught/entangled

### Integral monitoring line

Optimal positioning ensures consistent and accurate monitoring



### Push and Twist connectors

Ensures a secure, leak free fit every time



### Silver Knight™

Anti-microbial protection, reducing bacteria by 99.9%



## B Select your bagging limb if required



### Compact Bagging Limb

or

### Silver Knight Bagging Limb










Limbs are available in varying lengths from 0.8-3m, dependent on the type of tubing chosen.

**A** + **B** = **A complete system ready for use**  
 Prioritising patient safety, providing clinical choice

## Ordering information





Code	Description	Box Qty.
------	-------------	----------

### A Flextube™ Silver Knight™ Breathing Systems

2016016	22mm Flextube Silver Knight anti-microbial breathing system with 2L bag, and integral monitoring line, ≥ 1.6m	12	
2016024	22mm Flextube Silver Knight anti-microbial breathing system with 2L bag, and integral monitoring line, ≥ 2.4m	8	
2016032	22mm Flextube Silver Knight anti-microbial breathing system with 2L bag, and integral monitoring line, ≥ 3.2m	5	
2016036	22mm Flextube Silver Knight anti-microbial breathing system with 2L bag, and integral monitoring line, ≥ 3.6m	8	
2016040	22mm Flextube Silver Knight anti-microbial breathing system with 2L bag, and integral monitoring line, ≥ 4.0m	4	
2016048	22mm Flextube Silver Knight anti-microbial breathing system with 2L bag, and integral monitoring line, ≥ 4.8m	4	
2016056	22mm Flextube Silver Knight anti-microbial breathing system with 2L bag, and integral monitoring line, ≥ 5.6m	6	

Make an enquiry

### Uniflow™ Silver Knight™ Breathing Systems

2919016	30mm Uniflow Silver Knight anti-microbial coaxial breathing system with 2L bag, and integral monitoring line, ≥ 1.6m	8	
2919024	30mm Uniflow Silver Knight anti-microbial coaxial breathing system with 2L bag, and integral monitoring line, ≥ 2.4m	10	
2919032	30mm Uniflow Silver Knight anti-microbial coaxial breathing system with 2L bag, and integral monitoring line, ≥ 3.2m	15	
2919048	30mm Uniflow Silver Knight anti-microbial coaxial breathing system with 2L bag, and integral monitoring line, ≥ 4.8m	10	





Make an enquiry

### B Compact™ Bagging Limbs

2016115	22mm Compact bagging limb, ≥ 1.5m	100
2016120	22mm Compact bagging limb, ≥ 2.0m	100
2016130	22mm Compact bagging limb, ≥ 3.0m	100

Make an enquiry

### Flextube™ Silver Knight™ Bagging Limbs

2016208	22mm Flextube Silver Knight anti-microbial limb ≥ 0.8m	15	
2016212	22mm Flextube Silver Knight anti-microbial bagging limb, ≥ 1.2m	30	
2016216	22mm Flextube Silver Knight anti-microbial bagging limb, ≥ 1.6m	25	
2016224	22mm Flextube Silver Knight anti-microbial bagging limb, ≥ 2.4m	25	

Make an enquiry

 Silver Knight

IS5.37 • Issue 4 03.25



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK  
T: +44 (0)118 965 6300 info@intersurgical.com www.intersurgical.com



The manufacturer Intersurgical Ltd is certified to ISO 14001:2015, ISO 9001:2015, ISO 13485:2016 and MDSAP

**Please think before you print**  
Save energy and paper.  
If you must print this information sheet please print it double sided.