



# BiTrac™ MaxShield Select™

Paediatric total face mask for Non-invasive Ventilation



Critical Care ▪ Pulmodyne®



Quality, innovation and choice

Interact with us



[www.intersurgical.com.au](http://www.intersurgical.com.au)



## BiTrac™ MaxShield Select™

Our total face mask for paediatrics, the BiTrac MaxShield Select, avoids the nose and disperses pressure over a larger surface area while reducing claustrophobia for the patient. Our interchangeable elbows greatly reduce waste since you can easily remove and replace the mask's elbow when transitioning patients between standard ventilators and bilevel machines.



### 3 sizes available

small adult, extra small and extra extra small

### Super soft touch bonnet headgear

### Secure clips

### Soft silicone seal

### Bayonet style lock

for secure elbow connection

### Interchangeable elbows available

to fit your patient's needs


### XS mask with vented elbow with anti-asphyxiation valve and extender tube



### XXS mask with standard non-vented elbow




### BiTrac™ MaxShield Select™ paediatric total face mask with interchangeable vented elbow with anti-asphyxiation valve

Code	Description	Box Qty.
313-9527	BiTrac MaxShield Select paediatric total face mask with interchangeable vented elbow with anti-asphyxiation valve, 15 cm extender tube with 22M swivel connection, extra small	10
 313-9526	BiTrac MaxShield Select paediatric total face mask with interchangeable vented elbow with anti-asphyxiation valve, 15 cm extender tube with 22M swivel connection, extra extra small	10




### BiTrac MaxShield Select paediatric total face mask with interchangeable non-vented standard elbow

Code	Description	Box Qty.
 313-9525	BiTrac MaxShield Select paediatric total face mask with interchangeable standard non-vented elbow, extra small	10
313-9524	BiTrac MaxShield Select paediatric total face mask with interchangeable standard non-vented elbow, extra extra small	10





### BiTrac MaxShield Select total face mask with interchangeable elbow, small adult

Code	Description	Box Qty.
313-9554	BiTrac MaxShield Select total face mask with interchangeable vented elbow with anti-asphyxiation valve, small adult	10
 313-9551	BiTrac MaxShield Select total face mask with interchangeable standard non-vented elbow, small adult	10



### Accessories

Code	Description	Box Qty.
313-9529	Replacement headgear for paediatric BiTrac MaxShield, size extra small	10
313-9528	Replacement headgear for paediatric BiTrac MaxShield, size extra extra small	10
 313-9602	BiTrac Select paediatric interchangeable vented elbow with anti-asphyxiation valve and 15 cm extender tube with 22M swivel connection	10
313-9599	BiTrac Select interchangeable standard non-vented elbow	10



### Size guides included

 illustrated

BiTrac MaxShield Select is manufactured by Pulmodyne Inc, USA and distributed by Intersurgical Ltd

Please think before you print Save energy and paper. If you must print this information sheet please print it double sided.

IS1.61 AU ▪ Issue 1 03.22



Intersurgical Australia Pty Ltd, 4/151 Beauchamp Road, Matraville, NSW, 2036, Australia  
T: +61 (0)2 8048 3300 info@intersurgical.com.au www.intersurgical.com.au

UK ▪ Ireland ▪ France ▪ Germany ▪ Spain ▪ Portugal ▪ Italy ▪ Benelux ▪ Sweden ▪ Denmark ▪ Lithuania ▪ Russia ▪ Czech Republic  
Turkey ▪ South Africa ▪ China ▪ Japan ▪ Taiwan ▪ Philippines ▪ USA ▪ Canada ▪ Colombia ▪ Australia