




O2-MAX™

Disposable CPAP device




Critical Care ▪ Pulmodyne®

O2-MAX™

Continuous positive airway pressure is becoming the standard of care for patients suffering from respiratory distress. The O2-MAX System can deliver up to 140 L/min of flow to the patient, with 30% FiO₂. Customise your FiO₂ to increase from 30, 60, or 90% with the Trio adaptor. Adjustable PEEP

options includes a 3-SET™ (5/7.5/10 cmH₂O) or 5-SET™ (5/7.5/10/12.5/15 cmH₂O) to suit a variety of patient needs. Additionally, this completely disposable system offers integrated nebulisation while maintaining a consistent flow to the patient.



O2-MAX accessories

Code	Description	Box Qty.
313-8000	O2-MAX Trio control unit	10
313-8300	O2-MAX exhalation filter	10
313-7716	O2-CPAP™ fixed value PEEP valve, 5.0 cmH ₂ O, blue	10
313-7717	O2-CPAP fixed value PEEP valve, 7.5 cmH ₂ O, yellow	10
313-7718	O2-CPAP fixed value PEEP valve, 10.0 cmH ₂ O, green	10
313-7719	O2-CPAP fixed value PEEP valve, 12.5 cmH ₂ O, pink	10
313-7720	O2-CPAP fixed value PEEP valve, 15.0 cmH ₂ O, red	10
313-7721	O2-CPAP fixed value PEEP valve, 20.0 cmH ₂ O, gold	10
313-7725	3-SET adjustable PEEP valve (5/7.5/10 cmH ₂ O), white	10
313-7730	5-SET adjustable PEEP valve (5/7.5/10/12.5/15 cmH ₂ O), grey	10



O2-MAX™ Fixed System with expandable tube (delivers 30% FiO₂)

Code	Description	Box Qty.
The systems below are provided with 3-SET™ adjustable PEEP valve		
313-7555XN	O2-MAX Fixed System with expandable tube, 3-SET adjustable PEEP valve, nebuliser and BiTrac ED Mask, large adult	10
313-7556XN	O2-MAX Fixed System with expandable tube, 3-SET adjustable PEEP valve, nebuliser and BiTrac ED Mask, medium adult	10
313-7555X	O2-MAX Fixed System with expandable tube, 3-SET adjustable PEEP valve and BiTrac ED Mask, large adult	10
313-7556X	O2-MAX Fixed System with expandable tube, 3-SET adjustable PEEP valve and BiTrac ED Mask, medium adult	10
The systems below are provided with 5-SET™ adjustable PEEP valve		
313-7554XN	O2-MAX Fixed System with expandable tube, 5-SET adjustable PEEP valve, nebuliser and BiTrac ED Mask, large adult	10
313-7557XN	O2-MAX Fixed System with expandable tube, 5-SET adjustable PEEP valve, nebuliser and BiTrac ED Mask, medium adult	10
313-7554X	O2-MAX Fixed System with expandable tube, 5-SET adjustable PEEP valve and BiTrac ED Mask, large adult	10
313-7557X	O2-MAX Fixed System with expandable tube, 5-SET adjustable PEEP valve and BiTrac ED Mask, medium adult	10

O2-MAX Trio System (delivers adjustable 30%, 60%, 90% FiO₂)

Code	Description	Box Qty.
The systems below are provided with expandable tube and nebuliser		
313-8018XN	O2-MAX Trio System with expandable tube, 5-SET adjustable PEEP valve, nebuliser and BiTrac ED Mask, large adult	10
313-8008XN	O2-MAX Trio System with expandable tube, 5-SET adjustable PEEP valve, nebuliser and BiTrac ED Mask, medium adult	10
The systems below are provided with smoothbore tube		
313-8018	O2-MAX Trio System with smoothbore tube, 5-SET adjustable PEEP valve and BiTrac ED Mask, large adult	10
313-8008	O2-MAX Trio System with smoothbore tube, 5-SET adjustable PEEP valve and BiTrac ED Mask, medium adult	10

BiTrac ED Masks with standard and CAP headgear

Code	Description	Box Qty.
313-7031	BiTrac ED Mask, extra large adult	10
313-7030	BiTrac ED mask, large adult	10
313-7029	BiTrac ED mask, medium adult	10
313-7028	BiTrac ED mask, small adult	10
313-7031C	BiTrac ED mask with CAP headgear, extra large adult	10
313-7030C	BiTrac ED mask with CAP headgear, large adult	10
313-7029C	BiTrac ED mask with CAP headgear, medium adult	10
313-7028C	BiTrac ED mask with CAP headgear, small adult	10

Quik-Connect™ disposable oxygen connections



Code	Description	Box Qty.
313-111111	Quik-Connect Ohmeda gas fitting, green	30
313-121111	Quik-Connect Chemetron gas fitting, green	30
313-131111	Quik-Connect Puritan Bennett gas fitting, green	30
313-141111	Quik-Connect Schrader gas fitting, green	30
313-351111	Quik-Connect American gas fitting, white	30
313-361111	Quik-Connect DIN gas fitting, white	30
313-371111	Quik-Connect AFNOR gas fitting, white	30
313-381111	Quik-Connect BS gas fitting, white	30
313-391111	Quik-Connect UNI gas fitting, white	30

3-SET, 5-SET, BiTrac, O2-MAX and Quik-Connect are manufactured by Pulmodyne Inc, USA and distributed by Intersurgical Ltd
 Ohmeda is a registered trademark of GE Healthcare Finland Oy. Puritan Bennett is a registered trademark of Covidien LP.
 Chemetron is a registered trademark of Allied Healthcare Products, Inc.

Please think before you print Save energy and paper. If you must print this information sheet please print it double sided.

IS1.64 AU ▪ Issue 1 03.22



Intersurgical Australia Pty Ltd, 4/151 Beauchamp Road, Matraville, NSW, 2036, Australia
 T: +61 (0)2 8048 3300 info@intersurgical.com.au www.intersurgical.com.au

UK ▪ Ireland ▪ France ▪ Germany ▪ Spain ▪ Portugal ▪ Italy ▪ Benelux ▪ Sweden ▪ Denmark ▪ Lithuania ▪ Russia ▪ Czech Republic
 Turkey ▪ South Africa ▪ China ▪ Japan ▪ Taiwan ▪ Philippines ▪ USA ▪ Canada ▪ Colombia ▪ Australia