

New oxygen recovery T-piece with reservoir bag

Delivers supplementary oxygen to post operative, spontaneously breathing patients, in transport and recovery.



Quality, innovation and choice

Interact with us



www.intersurgical.co.uk

New oxygen recovery T-piece with reservoir bag

We have expanded our range of oxygen recovery devices with the introduction of the new Oxygen recovery T-piece with reservoir bag.

This new device delivers supplementary oxygen to post operative, spontaneously breathing patients, in transport and recovery whilst the airway is still in place. The lightweight design means there is minimal pull on the patients supraglottic airway or tracheal tube and the reservoir bag provides an audible and visual indication of respiration.

Contact us for a sample or further information:

info@intersurgical.co.uk

For further information on this product visit:

www.intersurgical.co.uk/products/oxygen-and-aerosol-therapy/oxygen-recovery-t-piece-with-reservoir-bag



Code	Description	Tube length	Box Qty.
1041000	Oxygen recovery T-piece with reservoir bag		100
1041001	Oxygen recovery T-piece with reservoir bag and tube	2.1m	40
1041002	Oxygen recovery T-piece with reservoir bag and Intersurgical EcoLite™, adult, medium concentration oxygen mask and tube	2.1m	55



Quality, innovation and choice

T: 0118 9656 300

info@intersurgical.co.uk

www.intersurgical.co.uk

