

Product Name	22mm Flextube™ breathing system, ≥ 1.6m sterile		
Reference number	2000000S		
Product Type	Flextube Breathing System		
Manufacturer	Intersurgical		
Materials	 PP LDPE HDPE		
Packaging			
Unit Packaging	Individually packaged		
Box quantity	15		
Pallet quantity	120		

Instructions for Use	Instructions for use with product	Yes
	Instructions for use with each box	No


Applicable Standards	BS EN ISO 13485: Medical devices - Quality management systems - Requirements for regulatory purposes
	BS EN ISO 9001: Quality Management Systems - Requirements
	BS EN ISO 15223-1: Medical Devices - Symbols to be Used with Medical Device Labels, Labelling, & Information to be Supplied
	BS EN ISO 5367: Anaesthetic and Respiratory Equipment. Breathing sets and connectors
	BS EN ISO 5356-1: Anaesthetic and respiratory equipment. Conical connectors. Cones and sockets

Sterile	No
Shelf life	5 years
Disposal Methods	Following use, the product must be disposed of in line with local healthcare hygiene and waste disposal protocols and regulations
Re-use	 Single use only

Product characteristics

Technical information

Length (m)	≥ 1.6	
Resistance to flow at 30 L/min (cmH ₂ O)	Inspiratory	0.2
	Expiratory	0.2
Resistance to flow at 60 L/min (cmH ₂ O)	Inspiratory	0.7
	Expiratory	0.7
Compliance at 60 cmH ₂ O (ml/cmH ₂ O)	1.7	
Weight (g)	202	
Connectors (Machine - Patient)	22F - 22M/15F	



**Material
Abbreviations**

ABS	Acrylonitrile Butadiene Styrene
AL	Aluminium
Cu	Copper
EPDM	Ethylene Propylene Terpolymer Rubber
EVA	Ethylene-vinyl acetate
HDPE	High Density Polyethylene
HIPS	High Impact Polystyrene
LDPE	Low Density Polyethylene
LSR	Liquid Silicone Rubber
MBS	Methacrylate Butadiene Styrene
NRL	Natural Rubber Latex
PA	Polyamide (Nylon)
PC	Polycarbonate
PE	Polyethylene
IR	Polyisoprene rubber
POM	Polyoxymethylene
POP	Ethylene based plastomer
PP	Polypropylene
PPTA	Poly (p-phenylene terephthalamide)
PS	Polystyrene
PSU	Polysulfone
PU	Polyurethane
PVC	Poly Vinyl Chloride
SB	Styrene Butadiene
Si	Silicone
SS	Stainless Steel
TPE	Thermoplastic Elastomer
SBC	Styrene Butadiene Copolymer
SEBS	Styrene-Ethylene-Butadiene-Styrene

